

CONNECTION CABLES



BIONEN S.A.S.

Via P. Petrocchi, 42/1
50127 Firenze – Italia
www.bionen.it

BIONEN
medical devices

CONNECTION CABLES

GENERAL INFORMATION

Intended use: Transfer of biological signals picked up by the electrode to the recording machine.

Instructions for use: Connect the cable to the electrode and the recording machine.

Cleaning instructions: Wash with cold water and soap with pH 5,5.

Disinfection: Disinfect using validated processes either Amocid, Lysoformin, Alhyndex, IndicinGG or equivalent products, at high concentration too, without immersing the device into the liquids. Use products according to directives and regulations of the country in which the device is used.

Sterilization: NO STERILE device. If necessary, sterilize with validated processes for ETO or Gamma Ray sterilization.

Lifetime: After 5 years from sale, is not guaranteed the proper functionality of the device.

Storage: Particular precautions for storing are not necessary but it's best to put the electrodes in its box after washing.

Packaging: plastic pouch.

Contraindications and interactions: none.

Declared toxicity: none. LATEX FREE

Transport and Waste disposal: Any particular transport condition and waste disposal is expected.

Regulation: Class I device.

Comply with the MDD 93/42/CEE (D.L. 46/97) and manufactured according to the procedures of BIONEN Quality System certified ISO 13485.



AVAILABLE CONFIGURATIONS:

REF.	UM	N°	USE WITH:	L	PATIENT SIDE	MACHINE SIDE	TYPE	COLOR
0004.005	NR.	1	BRIDGE ELECTRODE	100 cm	CROCODILE	MCH D=2mm	MULTIPOLAR	GREY
0004.006	CF.	24	BRIDGE ELECTRODE	100 cm	CROCODILE	MCH D=2mm	MULTIPOLAR	GREY
0004.007	NR.	1	BRIDGE ELECTRODE	150 cm	SPRING MCH D=3mm	FMM D=1,5mm	MULTIPOLAR	GREY
0004.008	NR.	23	BRIDGE ELECTRODE	90 cm	CROCODILE	FMM D=1,5mm	MULTIPOLAR	GREY
0004.009	NR.	1	BRIDGE ELECTRODE	100 cm	Economic CROCODILE	Economic FMM D=1,5mm	MULTIPOLAR	GREY
0004.009-K	CF.	3	BRIDGE ELECTRODE	20 cm	COCCODRILLO Economico	FMM D=1,5mm Economico	MONOPOLARE	MULTICOLORE
0004.010	CF.	24	BRIDGE ELECTRODE	100 cm	CROCODILE	FMM D=1,5mm	MULTIPOLAR	GREY
0004.011	NR.	1	BRIDGE ELECTRODE	100 cm	SPRING MCH D=2mm	FMM D=1,5mm	MULTIPOLAR	GREY
0004.011-E	NR.	1	BRIDGE ELECTRODE	100 cm	Economic SPRING MCH D=2mm	Economic FMM D=1,5mm	MULTIPOLAR	GREY
0004.012	NR.	1	BRIDGE ELECTRODE	100 cm	MCH D=2mm	FMM D=1,5mm	MULTIPOLAR	GREY
0004.013	NR.	1	BRIDGE ELECTRODE	100 cm	SPRING MCH D=3mm	MCH D=2mm	MULTIPOLAR	GREY
0004.014	NR.	1	BRIDGE ELECTRODE	100 cm	SPRING MCH D=2mm	MCH D=2mm	MULTIPOLAR	GREY
0004.016	NR.	1	BRIDGE ELECTRODE	100 cm	MCH D=1,5mm	FMM D=1,5mm	MULTIPOLAR	GREY
0004.016-Y	NR.	1	BRIDGE ELECTRODE	100 cm	2 MCH D=1,5mm	FMM D=1,5mm	MULTIPOLAR	GREY
0004.016-200	NR.	1	BRIDGE ELECTRODE	200 cm	MCH D=1,5mm	FMM D=1,5mm	MULTIPOLAR	GREY
0004.017	NR.	1	BRIDGE ELECTRODE	100 cm	FMM D=1,5mm	FMM D=1,5mm	MULTIPOLAR	GREY
0004.018	NR.	1	BRIDGE ELECTRODE	100 cm	MCH D=4mm	FMM D=1,5mm	MULTIPOLAR	GREY
0004.019	NR.	1	ADAPTATOR	15 cm	MCH D=1,5mm	2 PIN KEY	MULTIPOLAR	MIX
0004.020	NR.	1	ADAPTATOR	300 cm	MCH D=1,5mm	2 PIN KEY	MULTIPOLAR	MIX
0004.030	NR.	1	NEMUS - SCM	50 cm	HR 212-10P-8P	DF13-4S-1.25C	MULTIPOLAR	BLACK
0004.031	NR.	1	SCM - HAND GRIP	200 cm	JACK 2.5MM 4 POLI	DF13-6S-1.25C	MULTIPOLAR	BLACK
0004.042	NR.	1	HEADCAPS 0013.XXX	500 cm	36 PIN MDR	36 PIN MDR	MULTIPOLAR	WHITE
0004.050-NK	NR.	1	EEG NIHION KODEN	400 cm	5 MCH D=1,5mm	50 PIN	MULTIPOLAR	MULTICOLOR
0004.060	NR.	1	D.I.B. - BluNet	10 cm	-	3 PIN KEY	MULTIPOLAR	GREY
0004.061	NR.	1	INDUCTIVE BELTS	15 cm	2 PIN KEY	2 PIN KEY	MULTIPOLAR	BLACK
0004.062	NR.	1	INDUCTIVE BELTS	60 cm	2 PIN KEY	2 PIN KEY	MULTIPOLAR	BLACK
0004.100	NR.	1	ST MARK PUDENDO	200 cm	4 PIN	2 TUCKEL 8/5 PIN	MULTIPOLAR	MIX
0004.201	NR.	1	EBNEURO SD8	150 cm	7 PIN	7 PIN	MULTIPOLAR	MIX
0004.202	NR.	1	EBNEURO SD8 (CPAP)	250 cm	5 PIN	RJ MASK	MULTIPOLAR	GREY
0004.203	NR.	1	BE MICRO	500 cm	6 PIN	6 PIN	MULTIPOLAR	BLACK
0004.204	NR.	1	BE MICRO	250 cm	USB	4 PIN	MULTIPOLAR	BLACK
0004.205	NR.	1	BE MICRO NONIN CABLE	125 cm	5 PIN	6 PIN	MULTIPOLAR	GREY
0004.206	NR.	1	LINK CABLE 256CH	54 cm	6 PIN	6 PIN	MULTIPOLAR	BLACK
0004.207	NR.	1	LINK CABLE 128CH	22 cm	6 PIN	6 PIN	MULTIPOLAR	BLACK
0004.208	NR.	1	BE PLUS LTM/GWi	150 cm	HR 212-10P-8P	DF13-6S-1.25C	MULTIPOLAR	MIX
0004.209	NR.	1	OXIMETER GWi	150 cm	HR 212-10P-8P	DF13-4S-1.25C	MULTIPOLAR	MIX
0004.213	NR.	1	BE MICRO	300 cm	6 PIN	6 PIN	MULTIPOLAR	BLACK
0004.215	NR.	1	BE MICRO FLS INTERFACE	200 cm	DF3A-7P-2DSA	DF13-6S-1.25C	MULTIPOLAR	MIX
0004.300	NR.	1	DISPOSABLE SURFACE ELECTRODES	100 cm	CROCODILE	FMM D=1,5mm	MULTIPOLAR	MULTICOLOR
0004.301	NR.	1	TESTINA POLIGRAFICA MONOMETRIA	250 cm	MDR 36 PIN	MDR 36 PIN	MULTIPOLAR	BLACK

CONNECTION CABLES

REF.	UM	N°	USE WITH:	L	PATIENT SIDE	MACHINE SIDE	TYPE	COLOR
0004.302	NR.	1	PRESSURE TRANSDUCER	250 cm	MDR 20 PIN	DF13-20DP-1.25C	MULTIPOLAR	BLACK
0004.303	NR.	1	NEMUS 2 PROXY	100 cm	36 PIN VHDCI	DF13-20DP-1.25C	MULTIPOLAR	BLACK
0004.304	NR.	1	DSM	90 cm	68 PIN VHDCI	68 PIN VHDCI	MULTIPOLAR	MIX
0004.305	NR.	1	DSM	500 cm	68 PIN VHDCI	68 PIN VHDCI	MULTIPOLAR	MIX
0004.306	NR.	1	STIMULATOR DSM	300 cm	3 PIN KEY	2 FMM D=1,5mm	MULTIPOLAR	MIX
0004.307	NR.	1	STIMULATOR DSM	150 cm	3 PIN KEY	2 FMM D=1,5mm	MULTIPOLAR	MIX
0004.308	NR.	1	LI-POL BATTERY CHARGER	500 cm	PP-018	-	MULTIPOLAR	MIX
0010.002	NR.	1	-	20 cm	-	FMM D=1,5mm/ MCH D=1,5mm	PONTICELLO	BLACK
0010.003	NR.	1	-	20 cm	-	MCH D= 2mm	PONTICELLO	BLACK
0010.004	CF.	2	-	20 cm	-	MCH D= 2mm	PONTICELLO	BLACK
0010.005	CF.	5	-	20 cm	-	FMM D=1,5mm/ MCH D=1,5mm	PONTICELLO	BLACK
0015.001	NR.	1	VEGA HEADCAPS 10 ELECTRODES	40 cm	FMM D=1,5mm	25 PIN	MULTIPOLAR	MULTICOLOR
0015.004	NR.	1	VEGA HEADCAPS 10 ELECTRODES	40 cm	MCH D=2mm	25 PIN	MULTIPOLAR	MULTICOLOR
0015.011	NR.	1	UNIVERSAL HEADCAPS 10 ELECTRODES	40 cm	FMM D=1,5mm	25 PIN	MULTIPOLAR	MULTICOLOR
0015.014	NR.	1	UNIVERSAL HEADCAPS 10 ELECTRODES	40 cm	MCH D=2mm	25 PIN	MULTIPOLAR	MULTICOLOR
0016.004	NR.	1	VEGA HEADCAPS 21 ELECTRODES	40 cm	FMM D=1,5mm	25 PIN	MULTIPOLAR	MULTICOLOR
0016.007	NR.	1	VEGA HEADCAPS 21 ELECTRODES	40 cm	MCH D=2mm	25 PIN	MULTIPOLAR	MULTICOLOR
0016.009	NR.	1	MICROMED HEADCAPS 21 ELECTRODES	40 cm	FMM D=1,5mm	25 PIN	MULTIPOLAR	MULTICOLOR
0016.011	NR.	1	UNIVERSAL HEADCAPS 21 ELECTRODES	40 cm	FMM D=1,5mm	25 PIN	MULTIPOLAR	MULTICOLOR
0016.011-64	NR.	1	UNIVERSAL HEADCAPS 64 ELECTRODES	40 cm	FMM D=1,5mm	3 CONNETTORI A 25 PIN	MULTIPOLAR	MULTICOLOR
0016.012	NR.	1	UNIVERSAL HEADCAPS 21 ELECTRODES	40 cm	MCH D=2mm	25 PIN	MULTIPOLAR	MULTICOLOR
0016.013	NR.	1	UNIVERSAL HEADCAPS 32+2 ELECTRODES + TP9 E TP10	40 cm	FMM D=1,5mm	37 PIN	MULTIPOLAR	MULTICOLOR
0016.014	NR.	1	ECI HEADCAPS 21 ELECTRODES	40 cm	FMM D=1,5mm	25 PIN	MULTIPOLAR	MULTICOLOR
0016.017	NR.	1	ECI HEADCAPS 21 ELECTRODES	40 cm	MCH D=2mm	25 PIN	MULTIPOLAR	MULTICOLOR
0016.111	NR.	1	UNIVERSAL HEADCAPS 32+2 ELECTRODES	40 cm	FMM D=1,5mm	37 PIN	MULTIPOLAR	MULTICOLOR
0016.112	NR.	1	UNIVERSAL HEADCAPS 61+2 ELECTRODES	40 cm	FMM D=1,5mm	2 CONNETTORI A 37 PIN	MULTIPOLAR	MULTICOLOR
0025.004	NR.	1	DISCONNECTIBLE MOPOLAR NEEDLE	100 cm	MCH D=2mm	2 PIN KEY	MULTIPOLAR	BLACK
0025.005	NR.	1	DISCONNECTIBLE MOPOLAR NEEDLE	100 cm	MCH D=0.7mm	2 PIN KEY	MULTIPOLAR	BLACK
0025.006	NR.	1	DISCONNECTIBLE MOPOLAR NEEDLE	120 cm	FMM D=1mm	MCH D=2mm	MULTIPOLAR	MIX
0025.007	NR.	1	DISCONNECTIBLE MOPOLAR NEEDLE	120 cm	FMM D=1mm	FMM D=1,5mm	MULTIPOLAR	MIX
0025.011	NR.	1	DISCONNECTIBLE MOPOLAR NEEDLE	120 cm	FMM D=1mm	MCH D=2mm	MULTIPOLAR	MULTICOLOR
0025.012	NR.	1	DISCONNECTIBLE MOPOLAR NEEDLE	120 cm	FMM D=1mm	FMM D=1,5mm	MULTIPOLAR	MULTICOLOR
0025.021	NR.	1	DISCONNECTIBLE MOPOLAR NEEDLE	120 cm	FMM D=1mm	MCH D=2mm	MULTIPOLAR	MULTICOLOR
0025.022	NR.	1	DISCONNECTIBLE MOPOLAR NEEDLE	120 cm	FMM D=1mm	FMM D=1,5mm	MULTIPOLAR	MULTICOLOR
0025.030	NR.	1	DISCONNECTIBLE MOPOLAR NEEDLE	100 cm	FMM D=1mm	MCH D=2mm	MULTIPOLAR	MULTICOLOR
0025.031	NR.	1	DISCONNECTIBLE MOPOLAR NEEDLE	100 cm	FMM D=1mm	FMM D=1,5mm	MULTIPOLAR	MULTICOLOR
0025.032	NR.	1	DISCONNECTIBLE MOPOLAR NEEDLE	100 cm	FMM D=0,7mm	FMM D=1,5mm	MULTIPOLAR	BLACK
0025.033	NR.	1	DISCONNECTIBLE MOPOLAR NEEDLE	100 cm	FMM D=0,7mm	FMM D=1,5mm	MULTIPOLAR	BROWN
0025.034	NR.	1	DISCONNECTIBLE MOPOLAR NEEDLE	100 cm	FMM D=0,7mm	FMM D=1,5mm	MULTIPOLAR	RED
0025.035	NR.	1	DISCONNECTIBLE MOPOLAR NEEDLE	100 cm	FMM D=0,7mm	FMM D=1,5mm	MULTIPOLAR	ORANGE

CONNECTION CABLES

REF.	UM	N°	USE WITH:	L	PATIENT SIDE	MACHINE SIDE	TYPE	COLOR
0025.036	NR.	1	DISCONNECTIBLE MOPOLAR NEEDLE	100 cm	FMM D=0,7mm	FMM D=1,5mm	MULTIPOLAR	YELLOW
0025.037	NR.	1	DISCONNECTIBLE MOPOLAR NEEDLE	100 cm	FMM D=0,7mm	FMM D=1,5mm	MULTIPOLAR	GREEN
0025.038	NR.	1	DISCONNECTIBLE MOPOLAR NEEDLE	100 cm	FMM D=0,7mm	FMM D=1,5mm	MULTIPOLAR	BLUE
0025.039	NR.	1	DISCONNECTIBLE MOPOLAR NEEDLE	100 cm	FMM D=0,7mm	FMM D=1,5mm	MULTIPOLAR	VIOLA
0025.040	NR.	1	DISCONNECTIBLE MOPOLAR NEEDLE	100 cm	FMM D=0,7mm	FMM D=1,5mm	MULTIPOLAR	GREY
0025.051	NR.	1	DISCONNECTIBLE MOPOLAR NEEDLE	100 cm	FMM D=0,7mm	FMM D=1,5mm	MULTIPOLAR	WHITE
0025.130	NR.	1	INSULATED DISCONNECTIBLE MOPOLAR NEEDLE	100 cm	FMM D=0,7mm	FMM D=1,5mm	MULTIPOLAR	MIX
0025.131	NR.	1	DISCONNECTIBLE MOPOLAR NEEDLE	150 cm	FMM D=0,7mm	FMM D=1,5mm	MULTIPOLAR	MIX
0025.132	NR.	1	DISPOSABLE SURFACE ELECTRODES	150 cm	FMM D=0,7mm	FMM D=1,5mm	MULTIPOLAR	GREY
0025.133	NR.	1	DISPOSABLE SURFACE ELECTRODES	200 cm	2 FMM D=1mm	2 FMM D=1,5mm	MULTIPOLAR	GREY
0025.137	NR.	1	DISPOSABLE SURFACE ELECTRODES	100 cm	2 FMM D=1mm	TUCKEL 5 PIN 240°	MULTIPOLAR	GREY
0025.138	NR.	1	DISPOSABLE SURFACE ELECTRODES	100 cm	FMM D=0,7mm	FMM D=1,5mm	MULTIPOLAR	GREY
0025.139	NR.	1	DISPOSABLE SURFACE ELECTRODES	200 cm	FMM D=0,7mm	FMM D=1,5mm	MULTIPOLAR	GREY
0025.140	NR.	1	DISPOSABLE SURFACE ELECTRODES	150 cm	FMM D=0,7mm	TUCKEL 5 PIN	MULTIPOLAR	GREY
0025.201	NR.	1	DISPOSABLE SURFACE ELECTRODES	100 cm	FMM D=0,7mm	TUCKEL 5 PIN	MULTIPOLAR	MIX
0025.202	NR.	1	DISPOSABLE SURFACE ELECTRODES	100 cm	FMM D=0,7mm	TUCKEL 7 PIN	MULTIPOLAR	MIX
0025.205	NR.	1	DISCONNECTIBLE MOPOLAR NEEDLE	150 cm	BIPOlar	FMM D=1,5mm	MULTIPOLAR	MIX
0025.301	NR.	1	DISCONNECTIBLE MOPOLAR NEEDLE	50 cm	BIPOlar	MCH D=2mm	MULTIPOLAR	MIX
0025.400	NR.	1	DISCONNECTIBLE MOPOLAR NEEDLE	60 cm	BIPOlar	TUCKEL 7 PIN	MULTIPOLAR	GREY
0025.401	NR.	1	DISCONNECTIBLE MOPOLAR NEEDLE	150 cm	FMM D=1mm	FMM D=1,5mm	MULTIPOLAR	RED/BLACK
0025.402	NR.	1	DISCONNECTIBLE MOPOLAR NEEDLE	150 cm	FMM D=1mm	FMM 15 PIN	MULTIPOLAR	MIX
0025.403	NR.	1	DISCONNECTIBLE MOPOLAR NEEDLE	150 cm	FMM D=1,5mm	DB 15 PIN PLUG	MULTIPOLAR	MIX
0044.001	NR.	1	DISCONNECTIBLE MOPOLAR NEEDLE	120 cm	DB 15 PIN PLUG	VHDCI 36 PIN PLUG	MULTIPOLAR	GREY
0044.002	NR.	1	DISCONNECTIBLE MOPOLAR NEEDLE	120 cm	DB 25 PIN PLUG	VHDCI 36 PIN PLUG	MULTIPOLAR	GREY
0044.003	NR.	1	SNAP CONNECTION	120 cm	SNAP	FMM D=1,5mm	MULTIPOLAR	RED
0044.004	NR.	1	SNAP CONNECTION	120 cm	SNAP	FMM D=1,5mm	MULTIPOLAR	BLACK
0044.005	NR.	1	SNAP CONNECTION	120 cm	SNAP	FMM D=1,5mm	MULTIPOLAR	GREEN
0044.006	CF.	3	SNAP CONNECTION	120 cm	SNAP	FMM D=1,5mm	MULTIPOLAR	BLUE
0044.011	NR.	1	SNAP CONNECTION	120 cm	SNAP	FMM D=1,5mm	MULTIPOLAR	YELLOW
0044.012	NR.	1	SNAP CONNECTION	120 cm	SNAP	FMM D=1,5mm	MULTIPOLAR	MULTICOLOR
0044.013	NR.	1	SNAP CONNECTION	120 cm	SNAP	MCH D=2mm	MULTIPOLAR	RED
0044.014	NR.	1	SNAP CONNECTION	120 cm	SNAP	MCH D=2mm	MULTIPOLAR	BLACK
0044.015	NR.	1	SNAP CONNECTION	120 cm	SNAP	MCH D=2mm	MULTIPOLAR	GREEN
0044.016	CF.	3	SNAP CONNECTION	120 cm	SNAP	MCH D=2mm	MULTIPOLAR	BLUE
0044.030	NR.	1	SNAP CONNECTION	45 cm	SNAP	MCH D=2mm	MULTIPOLAR	YELLOW
0044.031	NR.	1	SNAP CONNECTION	120 cm	SNAP	MCH D=2mm	MULTIPOLAR	MULTICOLOR
0044.033	NR.	1	SNAP CONNECTION	250 cm	2 SNAP	2 PIN KEY	MULTIPOLAR	YELLOW-BLUE
0044.034	NR.	1	SNAP CONNECTION	50 cm	2 SNAP	2 PIN KEY	MULTIPOLAR	YELLOW-BLUE
0044.051	NR.	1	SNAP CONNECTION	60 cm	2 SNAP	2 PIN KEY	MULTIPOLAR	YELLOW-BLUE
0044.061	NR.	1	SNAP CONNECTION	100 cm	3 SNAP	3 PIN KEY	MULTIPOLAR	BLACK-RED-GREEN
0044.062	CF.	3	SNAP CONNECTION	100 cm	SNAP	FMM D=1,5mm	MULTIPOLAR	RED
0044.071	NR.	1	SNAP CONNECTION	150 cm	SNAP	FMM D=1,5mm	MULTIPOLAR	RED

CONNECTION CABLES

REF.	UM	N°	USE WITH:	L	PATIENT SIDE	MACHINE SIDE	TYPE	COLOR
0044.072	NR.	1	SNAP CONNECTION	150 cm	SNAP	FMM D=1,5mm	MULTIPOLAR	BLACK
0044.073	NR.	1	SNAP CONNECTION	150 cm	SNAP	FMM D=1,5mm	MULTIPOLAR	RED
0044.074	NR.	1	SNAP CONNECTION	150 cm	SNAP	FMM D=1,5mm	MULTIPOLAR	BLACK
0044.101	NR.	1	SNAP CONNECTION	200 cm	SNAP	FMM D=1,5mm	MULTIPOLAR	GREEN
0044.102	NR.	1	SNAP CONNECTION	200 cm	SNAP	FMM D=1,5mm	MULTIPOLAR	BLUE
0044.103	NR.	1	SNAP CONNECTION	200 cm	SNAP	FMM D=1,5mm	MULTIPOLAR	RED
0044.104	NR.	1	SNAP CONNECTION	200 cm	SNAP	FMM D=1,5mm	MULTIPOLAR	BLACK
0044.105	NR.	1	SNAP CONNECTION	200 cm	SNAP	FMM D=1,5mm	MULTIPOLAR	GREEN
0044.111	NR.	1	SNAP CONNECTION	200 cm	SNAP	FMM D=1,5mm	MULTIPOLAR	BLUE
0044.112	NR.	1	SNAP CONNECTION	200 cm	SNAP	FMM D=1,5mm	MULTIPOLAR	YELLOW
0044.113	NR.	1	SNAP CONNECTION	200 cm	SNAP	MCH D=2mm	MULTIPOLAR	RED
0044.114	NR.	1	SNAP CONNECTION	200 cm	SNAP	MCH D=2mm	MULTIPOLAR	BLACK
0044.115	NR.	1	SNAP CONNECTION	200 cm	SNAP	MCH D=2mm	MULTIPOLAR	GREEN
0044.201	NR.	1	SNAP CONNECTION	15 cm	SNAP	MCH D=2mm	MULTIPOLAR	BLUE
0044.301	NR.	1	SNAP CONNECTION	200 cm	SNAP	MCH D=2mm	MULTIPOLAR	YELLOW
0056.001	NR.	1	SNAP CONNECTION	150 cm	SNAP	DRITTO	MULTIPOLAR	RED
0056.003	NR.	1	GROUND ELECTRODE	150 cm	CLIP	FMM D=1,5mm	MULTIPOLAR	GREY
0056.004	NR.	1	ECG (GALILEO)	150 cm	SNAP	FMM D=1,5mm	MULTIPOLAR	GREY
0056.005	NR.	1	ECG (GALILEO)	150 cm	SNAP	TUCKEL 5 PIN	MULTIPOLAR	GREY
0056.006	NR.	1	ECG (GALILEO)	150 cm	SNAP	TUCKEL 5 PIN	MULTIPOLAR	GREY
0056.008	NR.	1	ECG (GALILEO)	150 cm	SNAP	FMM D=1,5mm	MULTIPOLAR	GREY
0057.001	NR.	1	ECG (GALILEO)	120 cm	SNAP	MCH D=2mm	MULTIPOLAR	GREY
0057.003	NR.	1	ECG (KEYPOINT)	120 cm	3 SNAP	TUCKEL 5 PIN	MULTIPOLAR	GREY
0057.004	NR.	1	EMG	120 cm	3 CROCODILE	TUCKEL 5 PIN	MULTIPOLAR	GREY
0057.011	NR.	1	EMG	120 cm	3 CROCODILE	3 FMM D=1,5mm	MULTIPOLAR	GREY
0057.012	NR.	1	EMG	120 cm	3 CROCODILE	3 MCH D=2mm	MULTIPOLAR	GREY
0057.103	NR.	1	EMG	120 cm	2 CROCODILE	TUCKEL 5 PIN	MULTIPOLAR	GREY
0057.104	NR.	1	EMG	100 cm	CROCODILE	TUCKEL 5 PIN	MULTIPOLAR	GREY
0057.106	NR.	1	EMG	120 cm	3 CROCODILE P1	TUCKEL 5 PIN	MULTIPOLAR	BLACK
0057.107	NR.	1	EMG	120 cm	2 CROCODILE	3 PIN KEY	MULTIPOLAR	BLACK
0058.001	NR.	1	EMG	350 cm	3 CROCODILE	3 FMM D=1,5mm	MULTIPOLAR	
0058.010	NR.	1	EEG/EMG	350 cm	-	OCCHIELLO - SPINA SHUCO	MULTIPOLAR	YELLOW / GREEN
0058.011	NR.	1	EEG/EMG	350 cm	-	OCCHIELLO - PINZA	MULTIPOLAR	YELLOW / GREEN
0058.012	NR.	1	EEG/EMG	350 cm	-	PINZA	MULTIPOLAR	YELLOW / GREEN
0058.021	NR.	1	EEG/EMG	500 cm	-	PINZA	MULTIPOLAR	YELLOW / GREEN
0058.022	NR.	1	EEG/EMG	350 cm	-	PINZA - CONNETTORE NODO EQUIPOTENZIALE	MULTIPOLAR	YELLOW / GREEN
0064.001	NR.	1	EEG/EMG	150 cm	-	PINZA - CONNETTORE NODO EQUIPOTENZIALE	MULTIPOLAR	YELLOW / GREEN
0064.002	NR.	1	ADAPTER	150 cm	TUCKEL 5 PIN FMM	TUCKEL 3 PIN MCH	BIPOLAR	GREY
0064.003	NR.	1	ADAPTER	4 cm	TUCKEL 5 PIN FMM	TUCKEL 3 PIN		
0064.010	NR.	1	ADAPTER	50 cm	TUCKEL 7 PIN MCH	TUCKEL 5 PIN FMM	BIPOLAR	GREY

CONNECTION CABLES

REF.	UM	N°	USE WITH:	L	PATIENT SIDE	MACHINE SIDE	TYPE	COLOR
0064.011	NR	1	NEMUS ADAPTER	30 cm	2 JACK MONO FMM 6,3 mm	4 PIN KEY	MULTIPOLAR	BLACK