

# SILICONE NET HEADCAP

BIONEN Silicon NET, with their design, provide stability to bridge electrodes during EEG examinations. Made with soft silicone, they ensure patient comfort and precise electrode placement. Available in various sizes (for adults and children) and with options of variable bands (4x2, 6x2, and 6x4), they can be equipped with chin straps in PELD or Neoprene.

## KEY Feature:

- High comfort for the patient
- Available with different number of bands (4x2, 6x2, 6x4).
- Available in adult and pediatric sizes.



## SPECIFICATIONS

SILICONE HEADCAP	
Sterilization	Not sterile
Shelflife	5 years (60 months)
Side Effects and interactions	None
Possibile toxicity declared	None. LATEX FREE

  

REGULATION	
CND Classification	Z12109099
Basic UDI	8033315120013RN
Risk class	CE - class I FDA – Approved by FDA through 510(k) - class I
Country of origin	Italy
Regulation conformity	MDR (UE) 2017/745 Attachment I
Built following procedures of SAQ BIONEN	ISO 13485:2016

NOTE: DM's cleaning, disinfection, manipulation and preservation modalities are included in Instructions for Use.

## MATERIALS

HEADCAP	
Headcap	Silicone
Chin strap	Riblene/Velcro
Intersections	ABS
Pipe stoppers	PVC
Tape	Cotton
Buttons	Nickel-plated brass 0/1 BOT

  

PACKAGING	
Box	1/box

## AVAILABLE CONFIGURATIONS - 0005-XXX

REF	QTY	N° ROWS	Ø SILICONE	CHIN STRAP	SIZE
0005.000	1	6x4	4mm	RIBLENE	ADULT
0005.001	1	6x2	4mm	RIBLENE	ADULT
0005.002	1	4x2	4mm	RIBLENE	PEDIATRIC
0005.003	1	6x2	4mm	VELCRO/LINER	ADULT
0005.004	1	4x2	4mm	VELCRO/LINER	PEDIATRIC
0005.005	1	4x2	3mm	WHITE VELCRO	PEDIATRIC
0005.006	1	6x4	4mm	VELCRO/LINER	ADULT/PEDIATRIC
0005.007	1	5x4	4mm	VELCRO/LINER	ADULT